































<b>Name of Pack: AO Modular Hand System</b>		<b>Specialty: Orthopaedics</b>	
<b>MaQs No: 72985</b>		<b>Distributed To: Sterile Store Room L2 MCC</b>	
<b>Number of Sets: 1</b> <b>DPW = Don't Pack Without.</b>		<b>Sterilisation Process: Steam</b>	
<b>Picture of Item</b>	<b>Quantity In Set</b>	<b>Name of Item</b>	<b>Cat. No.</b>
	<b>1</b> <b>DPW</b>	<b>Handle with mini quick coupling</b>	<b>311.012</b>
	<b>1</b> <b>DPW</b>	<b>1.5/1.1mm Double Drill Guide</b>	<b>312.140</b>
	<b>1</b> <b>DPW</b>	<b>2.4mm Screwdriver with Holding Sleeve</b>	<b>313.96</b> <b>313.97</b>
	<b>1</b> <b>DPW</b>	<b>Depth Gauge for 1.3/1.5mm Screws</b>	<b>319.004</b>
	<b>1</b> <b>DPW</b>	<b>Depth Gauge for 2.0/2.4mm Screws &amp; Protector</b>	<b>319.006</b>
	<b>1</b> <b>DPW</b>	<b>1.3/1.0mm Double Drill Sleeve</b>	<b>319.134</b>
	<b>1</b>	<b>Sharp Hook</b>	<b>319.390</b>
	<b>1</b> <b>DPW</b>	<b>2.0/1.5mm Universal Drill Guide</b>	<b>323.201</b>
	<b>1</b> <b>DPW</b>	<b>2.4/1.8mm Universal Drill Guide</b>	<b>323.202</b>
	<b>1</b>	<b>Screw &amp; Plate Holding Forceps</b>	<b>347.985</b>
	<b>1</b>	<b>Bending/Cutting Pliers</b>	<b>391.962</b>
	<b>2</b>	<b>Universal Bending Pliers</b>	<b>391.963</b>













<b>Sterile Services Unit – Waikato DHB</b>		
<b>Name of Set: AO Modular Hand System</b>		<b>Version 4.0</b>
<b>CNC Authorised:</b>	<b>SSU Authorised:</b>	<b>Date Reviewed: November 15</b>

	1	Reduction Forceps with Points wide Soft lock 127mm	399.074
	1	Plate Holding Forceps with Swivel Foot	398.819
	1	90mm Termite Forceps	398.950
	1	Small Hohman Retractor 6mm short narrow tip	399.180
	1	Small Hohman Retractor 8mm short narrow tip	399.190
	1	Periosteal Elevator	399.481
	1	Periosteal Elevator	399.480
	1	Reduction Forceps with points Soft lock 126mm	399.071
		Pack in own grey plastic trays and place in top of metal container	
1.0/1.3mm Hand Module Case		Self Tapping Titanium Cortex Screw	
	2	1.0 x 6mm T2 D11 S1	400.526
	2	1.0 x 7mm T2 D11 S2	400.527
	2	1.0 x 8mm T2 D11 S3	400.528
	2	1.0 x 9mm T2 D11 S4	400.529

Sterile Services Unit – Waikato DHB		
Name of Set: AO Modular Hand System		Version 4.0
CNC Authorised:	SSU Authorised:	Date Reviewed: November 15













	2	1.0 x 10mm <b>T2 D11 S5</b>	400.530
	2	1.0 x 11mm <b>T2 D11 S6</b>	400.531
	2	1.0 x 12mm <b>T2 D11 S7</b>	400.532
	2	1.0 x 13mm <b>T2 D11 S8</b>	400.533
	2	1.0 x 14mm <b>T2 D11 S9</b>	400.534
		<b>Self Tapping Titanium Cortex Screw</b>	
	4	1.3 x 6mm <b>T2 D11 S10</b>	400.686
	4	1.3 x 7mm <b>T2 D11 S11</b>	400.687
	4	1.3 x 8mm <b>T2 D11 S12</b>	400.688
	4	1.3 x 9mm <b>T2 D11 S13</b>	400.689
	4	1.3 x 10mm <b>T2 D11 S14</b>	400.690
	4	1.3 x 11mm <b>T2 D11 S15</b>	400.691
	4	1.3 x 12mm <b>T2 D11 S16</b>	400.692
	4	1.3 x 13mm <b>T2 D11 S17</b>	400.693
	4	1.3 x 14mm <b>T2 D11 S18</b>	400.694
	4	1.3 x 16mm <b>T2 D11 S19</b>	400.696
	4	1.3 x 18mm <b>T2 D11 S20</b>	400.698
	1	1.3mm Titanium Straight Plate 6 hole (23mm)	<b>T2 D12 S30</b> 421.306
	1	1.3mm Titanium Straight Plate 12 hole (47mm)	<b>T2 D12 S31</b> 421.312
	1	1.3mm Titanium Extended H-Plate Right (15mm)	<b>T2 D12 S32</b> 421.320
	1	1.3mm Titanium Extended H-Plate Left (15mm)	<b>T2 D12 S33</b> 421.321
	1	1.3mm Titanium T-Plate 3 hole x 8 hole (32mm)	<b>T2 D12 S34</b> 421.333
	1	1.3mm Titanium T-Plate 4 hole x 8 hole (32mm)	<b>T2 D12 S35</b> 421.334
	1	1.3mm Titanium Y-Plate 3 hole x 8 hole (32mm)	<b>T2 D12 S36</b> 421.335

Sterile Services Unit – Waikato DHB		
Name of Set: AO Modular Hand System		Version 4.0
CNC Authorised:	SSU Authorised:	Date Reviewed: November 15









	1	Countersink for 1.3/1.5mm Screws	310.971
	1 DPW	1.3mm Cruciform Screwdriver Blade with H/Sleeve	314.411
	1 DPW	1.0mm Cruciform Screwdriver Blade with H/Sleeve	314.480
	1 DPW	0.76mm Drill Bit with 6mm stop Mini Q.C.	316.286
	1	0.76mm Drill Bit with 8mm stop Mini Q.C.	316.288
	1	0.76mm Drill Bit with 10mm Stop Mini Q.C.	316.290
	1 DPW	0.76mm Drill Bit with 12mm Stop Mini Q.C.	316.292
	1	0.76 Drill Bit with 14mm stop Mini Q.C.	316.294
	2 DPW At least 1.	1.0mm Drill Bit Mini Q.C.	316.396
	2 DPW At least 1.	1.3mm Drill Bit Mini Q.C.	316.402
		<b>Place in screw caddy and place caddy in the side of metal container</b>	
1.5 Hand Module Case		Self Tapping Titanium Cortex Screw	
	4	1.5 x 6mm <b>T2 D11 S21</b>	400.806
	4	1.5 x 7mm <b>T2 D11 S22</b>	400.807
	4	1.5 x 8mm <b>T2 D11 S23</b>	400.808
	4	1.5 x 9mm <b>T2 D11 S24</b>	400.809
	4	1.5 x 10mm <b>T2 D11 S25</b>	400.810
	4	1.5 x 11mm <b>T2 D11 S26</b>	400.811
	4	1.5 x 12mm <b>T2 D11 S27</b>	400.812

Sterile Services Unit – Waikato DHB		
Name of Set: AO Modular Hand System		Version 4.0
CNC Authorised:	SSU Authorised:	Date Reviewed: November 15














	4	1.5 x 13mm <b>T2 D11 S28</b>	400.813
	4	1.5 x 14mm <b>T2 D11 S29</b>	400.814
	4	1.5 x 16mm <b>T2 D11 S30</b>	400.816
	4	1.5 x 18mm <b>T2 D11 S31</b>	400.818
	4	1.5 x 20mm <b>T2 D11 S32</b>	400.820
	4	1.5 x 22mm <b>T2 D11 S33</b>	400.822
	4	1.5 x 24mm <b>T2 D11 S34</b>	400.824
	1	1.5mm Titanium Straight Plate 6 hole (29mm)	446.031 <b>T2 D12 S39</b>
	1	1.5mm Titanium Straight Plate 12 hole (59mm)	446.032 <b>T2 D12 S40</b>
	1	1.5mm Titanium T-Plate 3 hole x 8 hole (40mm)	446.233 <b>T2 D12 S41</b>
	1	1.5mm Titanium T-Plate 4 hole x 8 hole (40mm)	446.234 <b>T2 D12 S42</b>
	1	1.5mm Titanium Extended H-Plate Right (18mm)	446.482 <b>T2 D12 S37</b>
	1	1.5mm Titanium Extended H-Plate Left (18mm)	446.483 <b>T2 D13 S1</b>
	1	1.5mm Titanium Y-Plate 3 hole x 1 hole	446.612 <b>T2 D13 S2</b>
	1	1.5mm Titanium Condylar Plate 7 hole Left (36mm)	446.630 <b>T2 D13 S3</b>
	1	1.5mm Titanium Condylar Plate 7 hole Right (36mm)	446.640 <b>T2 D13 S4</b>
	2 <b>DPW</b> At least 1.	1.1mm Drill Bit	310.111 <b>T2 D2 S16</b>
	2 <b>DPW</b> At least 1.	1.5mm Drill Bit	310.141 <b>T2 D2 S17</b>
	1	Countersink for 1.3/1.5mm Screws	310.971

Sterile Services Unit – Waikato DHB		
Name of Set: AO Modular Hand System		Version 4.0
CNC Authorised:	SSU Authorised:	Date Reviewed: November 15











	<b>1</b> <b>DPW</b>	<b>1.5/2.0mm Cruciform Screwdriver Blade with H/Sleeve</b>	<b>314.667</b>
		<b>Place in screw caddy and place caddy in the side of metal container</b>	
<b>2.0 Hand Module Case</b>  		<b>Self Tapping Titanium Cortex Screw</b>	
	<b>4</b>	<b>2.0 x 6mm T2 D11 S35</b>	<b>401.806</b>
	<b>4</b>	<b>2.0 x 7mm T2 D11 S36</b>	<b>401.807</b>
	<b>4</b>	<b>2.0 x 8mm T2 D11 S37</b>	<b>401.808</b>
	<b>4</b>	<b>2.0 x 9mm T2 D11 S38</b>	<b>401.809</b>
	<b>4</b>	<b>2.0 x 10mm T2 D11 S39</b>	<b>401.810</b>
	<b>4</b>	<b>2.0 x 11mm T2 D11 S40</b>	<b>401.811</b>
	<b>4</b>	<b>2.0 x 12mm T2 D11 S41</b>	<b>401.812</b>
	<b>4</b>	<b>2.0 x 13mm T2 D11 S42</b>	<b>401.813</b>
		<b>2.0 x 14mm T2 D12 S1</b>	<b>401.814</b>
	<b>4</b>	<b>2.0 x 16mm T2 D12 S2</b>	<b>401.816</b>
	<b>4</b>	<b>2.0 x 18mm T2 D12 S3</b>	<b>401.818</b>
	<b>4</b>	<b>2.0 x 20mm T2 D12 S4</b>	<b>401.820</b>
	<b>4</b>	<b>2.0 x 22mm T2 D12 S5</b>	<b>401.822</b>
	<b>4</b>	<b>2.0 x 24mm T2 D12 S6</b>	<b>401.824</b>
	<b>4</b>	<b>2.0 x 26mm T2 D12 S7</b>	<b>401.826</b>
	<b>2</b>	<b>2.0 x 28mm T2 D12 S8</b>	<b>401.828</b>
	<b>2</b>	<b>2.0 x 30mm T2 D12 S9</b>	<b>401.830</b>
	<b>2</b>	<b>2.0 x 32mm T2 D12 S10</b>	<b>401.832</b>
	<b>1</b>	<b>2.0mm Titanium Condylar Plate 7 hole left (39mm)</b>	<b>TIS</b> <b>443.610</b>
	<b>1</b>	<b>2.0mm Titanium Condylar Plate 7 hole Right (39mm)</b>	<b>443.620</b> <b>T2 D12</b> <b>S38</b>
	<b>1</b>	<b>2.0mm Titanium Straight Plate 6 hole (35mm)</b>	<b>447.031</b> <b>T2 D13 S5</b>
	<b>1</b>	<b>2.0mm Titanium Straight Plate 12 hole (71mm)</b>	<b>447.032</b> <b>T2 D13 S6</b>
	<b>1</b>	<b>2.0mm Titanium T-Plate 2 hole x 8 hole (48mm)</b>	<b>447.232</b> <b>T2 D13 S7</b>

Sterile Services Unit – Waikato DHB		
Name of Set: AO Modular Hand System		Version 4.0
CNC Authorised:	SSU Authorised:	Date Reviewed: November 15

	1	2.0mm Titanium T-Plate 3 hole x 8 hole (48mm)	447.233 <b>T2 D13 S8</b>
	1	2.0mm Titanium Y-Plate 3 hole x 8 hole (48mm)	<b>T2 D13 S9</b> 447.612
	1	2.0mm Titanium LC-DCP 4 hole (24mm)	<b>T2 D13 S21</b> 449.931
	1	2.0mm Titanium LC-DCP 6 hole (36mm)	<b>T2 D13 S22</b> 449.936
	1	2.0mm Titanium LC-DCP 8 hole (48mm)	<b>T2 D13 S23</b> 449.938
	2 DPW At least 1.	1.5mm Drill Bit	<b>T2 D2 S17</b> 310.141
	2 DPW At least 1.	2.0mm Drill Bit	310.201 or 513.140 <b>T1 D2</b>
	1	Countersink for 2.0/2.4mm screws	310.972
	1 DPW	1.5/2.0mm Cruciform Screwdriver Blade with H/sleeve	314.670
		<b>Place in screw caddy and place caddy in the side of metal container</b>	
2.4 Hand Module Case  	4 4 4 4 4 4 4 4 4	Self-Tapping Titanium Cortex Screw 2.4 x 6mm <b>T2 D12 S11</b> 2.4 x 7mm <b>T2 D12 S12</b> 2.4 x 8mm <b>T2 D12 S13</b> 2.4 x 9mm <b>T2 D12 S14</b> 2.4 x 10mm <b>T2 D12 S15</b> 2.4 x 11mm <b>T2 D12 S16</b> 2.4 x 12mm <b>T2 D12 S17</b> 2.4 x 13mm <b>T2 D12 S18</b> 2.4 x 14mm <b>T2 D12 S19</b>	401.506 401.507 401.508 401.509 401.510 401.511 401.512 401.513 401.514








Sterile Services Unit – Waikato DHB		
Name of Set: AO Modular Hand System		Version 4.0
CNC Authorised:	SSU Authorised:	Date Reviewed: November 15



	4	2.4 x 16mm <b>T2 D12 S20</b>	401.516
	4	2.4 x 18mm <b>T2 D12 S21</b>	401.518
	4	2.4 x 20mm <b>T2 D12 S22</b>	401.520
	4	2.4 x 22mm <b>T2 D12 S23</b>	401.522
	4	2.4 x 24mm <b>T2 D12 S24</b>	401.524
	2	2.4 x 26mm <b>T2 D12 S25</b>	401.526
	2	2.4 x 28mm <b>T2 D12 S26</b>	401.528
	2	2.4 x 30mm <b>T2 D12 S27</b>	401.530
	2	2.4 x 32mm <b>T2 D12 S28</b>	401.532
	4	1.8mm Titanium Buttress Pin (18mm)	401.968
	1	2.4mm Titanium Straight Plate 6 hole (36mm)	<b>T2 D13 S10</b> 449.906
	1	2.4mm Titanium Straight Plate 12 hole (72mm)	<b>T2 D13 S11</b> 449.912
	1	2.4mm Titanium T-Plate 2 hole x 8 hole (54mm)	<b>T2 D13 S13</b> 449.913
	1	2.4mm Titanium T-Plate 3 hole x 8 hole (54mm)	<b>T2 D13 S14</b> 449.914
	1	2.4mm Titanium Y-Plate 3 hole x 8 hole (56mm)	<b>T2 D13 S15</b> 449.915
	1	2.4mm Titanium Condylar Plate 8 hole Right (57mm)	<b>T2 D13 S16</b> 449.916
	1	2.4mm Titanium Condylar Plate 8 hole Left (57mm)	<b>T2 D13 S17</b> 449.917
	1	2.4mm Titanium LC-DCP 4 hole (31mm)	<b>T2 D13 S18</b> 449.924
	1	2.4mm Titanium LC-DCP 6 hole (47mm)	<b>T2 D13 S19</b> 449.926

Sterile Services Unit – Waikato DHB		
Name of Set: AO Modular Hand System		Version 4.0
CNC Authorised:	SSU Authorised:	Date Reviewed: November 15



	1	2.4mm Titanium LC-DCP 8 hole (63mm)	<b>T2 D13 S20 449.928</b>
	1	Countersink for 2.0/2.4mm Screws	310.972
	1	Threaded Drill Guide for 2.4mm Condylar Plate	312.920
	2 DPW At least 1.	1.8mm Drill Bit Mini Q.C.	317.861 <b>T2 D2 S18</b>
	2 DPW At least 1.	2.4mm Drill Bit Mini Q.C.	317.871
	1	Bending Iron for 2.4mm Condylar Plate	329.090
		Place in screw caddy and place caddy in the side of metal container	
<b>Packing Instructions</b>		<b>Pack in it's own metal container and lid. Use large tray liner. Use Pink Sterilisation Corner Protectors. Wrap in 121cm x 121cm Heavy Duty Quick Check.</b>	

Sterile Services Unit – Waikato DHB		
Name of Set: AO Modular Hand System		Version 4.0
CNC Authorised:	SSU Authorised:	Date Reviewed: November 15